

## Background and Objective

### Background

EMT (Epithelial-Mesenchymal transition) is a significant event in tumor metastasis and malignancy. Cancer cells are stimulated by EMT-inducible agonists from the surrounding stromal cells, and cause transformation such as weakening of cell-cell adhesion and the accompanying acceleration of cell migration and invasion. Inhibition of EMT is considered to enable controlling of malignant transformation.

### Study Objective

The objective of this study was development of a high-throughput assay system of EMT inhibitor screen in 3D cell culture.

## Methods

### 3D cell culture

We adopted 3D cell culture method on NanoCulture<sup>®</sup> Plate (SCIVAX), a scaffold type spheroid culture plate, as a culture method which can confirm spheroid morphologies.

### EMT induction and inhibition

TGF- $\beta$ 2 and its inhibitor, SB431542, were used as an EMT inducer and as a positive control for EMT inhibitor, respectively.

### Gene expression analysis

Expressions of E-cadherin, N-cadherin, vimentin, and zeb-1, represent EMT maker genes, were evaluated by qRT-PCR. TBP gene expressions was also evaluated as internal control.


### Intra-spheroid hypoxia sensing

The intra-spheroid hypoxic area was visualized using Hypoxia Probe (SCIVAX).

### Evaluation of the EMT status


Fluorescence intensity of the Hypoxia Probe of intra-spheroid in whole well was detected using Celigo<sup>®</sup> Imaging Cytometer (Brooks Automation, Inc.), as an indicator for EMT induction and inhibition.

## Results




**Figure 1. TGF- $\beta$ 2 treatment induces EMT, which is characterized by decreasing cell-cell adherence, increasing cell migration, and change of expression of EMT maker genes.**

(A) A549 spheroids treated with TGF- $\beta$ 2 (5 ng/mL) at day 3. Time-Lapse analysis revealed migration of cells from spheroids, leads to spheroid collapse. (B) SB431542 inhibited the spheroid collapse and change of expression of EMT maker genes.



**Figure 2. Fluorescence intensity of Hypoxia Probe reflects degree of cell-cell adhesion and spheroid morphology.**

(A) Comparison of E-cadherin expression of epithelial cancer cell lines. (B) Cell lines form spheroids of various sizes and shapes, according to these E-cadherin expression level. Differential changes of signal intensity of Hypoxia Probe depends on the spheroid features.



**Figure 3. Degree of intra-spheroid hypoxia reflects the spheroid compactness (EMT status).**

(A) Spheroids was treated with TGF- $\beta$ 2, or cotreated with TGF- $\beta$ 2 (5 ng/mL) and SB431542, for 3 days. Fluorescence intensity of Hypoxia Probe decreased according to TGF- $\beta$ 2 concentration, or increased according to SB431542 concentration (b). In these condition, cell viability was not affected. (B) Fluorescence intensity of Hypoxia Probe depended on the morphology of EMT-induced spheroids. The signal intensity was measured using Celigo<sup>®</sup> Imaging Cytometer.

## Conclusion

- We developed unique high-throughput assay system for EMT inhibitor screen, by evaluating intra-spheroid hypoxicity as an indicator for alternation of the spheroid morphology.
- We are now running a screen for EMT inhibitor candidate compounds, by using this assay system.